

- EN** Instructions for use
- FR** Mode d'emploi
- ES** Instrucciones de uso
- RU** Инструкция по применению
- RO** Instrucțiuni de utilizare
- EL** Οδηγίες χρήσης
- DE** Gebrauchsanleitung
- NL** Gebruiksaanwijzing
- IT** Istruzioni per l'uso
- PL** Instrukcja obsługi
- SV** Bruksanvisning
- PT** Instruções de uso

EC REP


SHOFU DENTAL GmbH
Am Brüll 17, 40878 Ratingen,
Germany

SHOFU DENTAL CORPORATION
1225 Stone Drive, San Marcos,
California 92078-4059, USA


SHOFU DENTAL ASIA-PACIFIC PTE.LTD.
10 Science Park Road, #03-12,
The Alpha, Science Park II, Singapore 117684



○ 2-layer technique



*** 3-layer technique**



COLOR CHART

■ Basic Shades

Shade	A1	A2	A3	A3.5	A4	rootA	B1	B2	B3	B4	rootB	C1	C2	C3	C4	rootC	D2	D3	D4
Opaque	A ₁ O	A ₂ O	A ₃ O	A _{3.5} O	A ₄ O	rootAO	B ₁ O	B ₂ O	B ₃ O	B ₄ O	rootBO	C ₁ O	C ₂ O	C ₃ O	C ₄ O	rootCO	D ₂ O	D ₃ O	D ₄ O
Opaque Dentin	OD-A ₁	OD-A ₂	OD-A ₃	OD-A _{3.5}	OD-A ₄	OD-rootA	OD-B ₁	OD-B ₂	OD-B ₃	OD-B ₄	OD-rootB	OD-C ₁	OD-C ₂	OD-C ₃	OD-C ₄	OD-rootC	OD-D ₂	OD-D ₃	OD-D ₄
Body	A ₁ B	A ₂ B	A ₃ B	A _{3.5} B	A ₄ B	rootAB	B ₁ B	B ₂ B	B ₃ B	B ₄ B	rootBB	C ₁ B	C ₂ B	C ₃ B	C ₄ B	rootCB	D ₂ B	D ₃ B	D ₄ B
Opal Porcelain Incisal	57	○					○*												
	58	*	○*					○*				○*	○*					○	
	59			○*	○*				○*					○*				*	○
	60					○*	○*			○*	○*				○*	○*		*	*
Translucent	T	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*

■ Red Shift Shades

Shade	R1	R2	R3	R3.5	R4	VR1	VR2	VR3	VR3.5	VR4
Opaque	R ₁ O	R ₂ O	R ₃ O	R _{3.5} O	R ₄ O	R ₁ O	R ₂ O	R ₃ O	R _{3.5} O	R ₄ O
Body	R ₁ B	R ₂ B	R ₃ B	R _{3.5} B	R ₄ B	VR ₁ B	VR ₂ B	VR ₃ B	VR _{3.5} B	VR ₄ B
Opal Porcelain Incisal	57	○				○	○			
	58	*	○*			*	*	○	○	
	59			○*	○*			*	*	○
	60				○*	○*				*
Translucent	T	*	*	*	*	*	*	*	*	*

■ Value Plus Shades

Shade	VA1	VA2	VA3	VA3.5	VA4	VB1	VB2	VB3	VB4
Opaque	A ₁ O	A ₂ O	A ₃ O	A _{3.5} O	A ₄ O	B ₁ O	B ₂ O	B ₃ O	B ₄ O
Body	VA ₁ B	VA ₂ B	VA ₃ B	VA _{3.5} B	VA ₄ B	VB ₁ B	VB ₂ B	VB ₃ B	VB ₄ B
Opal Porcelain Incisal	57	○	○			○	○		
	58	*	*	○	○	*	*	○	
	59			*	*	○		*	○
	60					*			*
Translucent	T	*	*	*	*	*	*	*	

■ Value Minus Shades

Shade	VmA1	VmA2	VmA3	VmA3.5	VmA4	VmrootA	VmR1	VmR2	VmR3	VmR3.5	VmR4	VmrootR
Opaque	VmA ₁ O	VmA ₂ O	VmA ₃ O	VmA _{3.5} O	VmA ₄ O	VmrootAO	VmR ₁ O	VmR ₂ O	VmR ₃ O	VmR _{3.5} O	VmR ₄ O	VmrootRO
Opaque Dentin	VmOD-A ₁	VmOD-A ₂	VmOD-A ₃	VmOD-A _{3.5}	VmOD-A ₄	VmOD-rootA	VmOD-R ₁	VmOD-R ₂	VmOD-R ₃	VmOD-R _{3.5}	VmOD-R ₄	VmOD-rootR
Body	VmA ₁ B	VmA ₂ B	VmA ₃ B	VmA _{3.5} B	VmA ₄ B	VmrootAB	VmR ₁ B	VmR ₂ B	VmR ₃ B	VmR _{3.5} B	VmR ₄ B	VmrootRB
Opal Porcelain Incisal	57	○					○					
	58	*	○*				*	○*				
	59			○*	○*				○*	○*		
	60					○*	○*				○*	○*
Translucent	T	*	*	*	*	*	*	*	*	*	*	

■ Whitening Shades

Shade	W1	W2	W3
Opaque	W ₁ O	W ₂ O	W ₃ O
Body	W ₁ B	W ₂ B	W ₃ B
Opal Porcelain Incisal	56	○	○ ^(Note 1)
	57	*	○ ^(Note 1) *
	58		
	59		
	60		
Translucent	*	*	*

Note 1: Mixing Ratios 56:57 → 2:1
Note 2: Mixing Ratios 56:57 → 1:2

FIRING SCHEDULE

Type of porcelain	Drying time	Temperature raising speed	Firing schedule
Opaque First firing Second firing	1-3 min.	50-60 °C/min	650 °C → 940-950 °C 0-1 min.
Margin First firing Second firing	3-5 min.	50-60 °C/min	650 °C → 940-950 °C 0-1 min.
Opaque Dentin, Body, Opal Porcelain, Incisal, Translucent First firing	5-7 min.	50-60 °C/min	650 °C → 910-930 °C 0-0.5 min.
Opaque Dentin, Body, Opal Porcelain, Incisal, Translucent Second firing	5-7 min.	50-60 °C/min	650 °C → 900-920 °C 0-0.5 min.
Add-on First firing	5-7 min.	50-60 °C/min	650 °C → 850-870 °C
Add-on (Glazing)	5-7 min.	50-60 °C/min	650 °C → 850-870 °C
Self-glazing	5-7 min.	50-60 °C/min	650 °C → 900-920 °C 0-0.5 min.

PHYSICAL PROPERTIES

	Type of firing	Coefficient of thermal expansion (25-500 °C)	Glass transition point
Opaque	2 nd firing	13.0×10 ⁻⁶ K ⁻¹	580 °C
	4 th firing		
Opaque Dentin	2 nd firing	13.0×10 ⁻⁶ K ⁻¹	580 °C
	4 th firing		
Body	2 nd firing	13.0×10 ⁻⁶ K ⁻¹	580 °C
	4 th firing		
Opal Porcelain Incisal	2 nd firing	13.0×10 ⁻⁶ K ⁻¹	580 °C
	4 th firing		
Translucent	2 nd firing	13.0×10 ⁻⁶ K ⁻¹	580 °C
	4 th firing		
Margin	2 nd firing	13.0×10 ⁻⁶ K ⁻¹	575 °C
	4 th firing		
Correction	2 nd firing	11.0×10 ⁻⁶ K ⁻¹	585 °C
	4 th firing		

: Vacuum firing (1.3-8.0 kPa)